# KN95 Respirator

# **Product Specifications**

### **Product Description**

KN95 Respirator with ear loops.

### **Basic Specifications**

- Item Number: HMA203WHJ2H
- Material: Melt Blown Cloth and Non-Woven Material
- Color: White
- Pack Size: 20

### **Additional Features**

- Disposable
- Non-Sterile
- 5-Ply
- Outside Layers: 2-Ply Non-Woven Material
- Middle Layer: Melt Blown Cloth
- Inner Layer: Latex Cotton
- Length of Earloop: 8"
- Size: Universal



**FOH** Health

FSSFNTIALS

#### **Quality Standard**

Click certification name to view

- FDA listed under product code MSH-Surgical Mask
- CE compliant under European Standard EN149:2001+A1:2009
- ASTM Level 1. BFE>=95%. Conforms to KN95, EN149; 2001 FFP2 standards.

#### Instructions

With sanitized hands, unfold the mask by pulling the sides of the mask apart. Put the mask over your nose and mouth. Adjust the straps to fit comfortably around the ears. Mold the metal bridge around your nose without pinching to make it fit snugly. Perform a seal check before use: cover the mask with both hands without pinching, and inhale and exhale. If you feel air blowing at the edge of the mask, readjust it. Correct fit of the mask requires contact with smooth skin. Discard after use.